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Literatur
1. Carraro P et al. Errors in a stat laboratory: Types and frequencies 10 years later. Clin Chem 2007; 53:7: 1338-42. 2. Joint Commission on Accreditation of Healthcare Organization. National Patient Goals. www.jointcommission.org Zuletzt abgerufen im Juni 2013. 3. Kohn S. Specimen mislabeling: A significant and costly cause of potentially serious medical errors. www.acutecaretesting.org April 2005. 4. CLSI. Procedures for the Collection of Arterial Blood Specimens; Approved Standard-Fourth Edition. CLSI document GP43-A4. Clinical and Laboratory Standards Institute, 950 West Valley Road, Suite 2500, Wayne, Pennsylvania 19087 USA, 2004. 5. CLSI. Blood Gas and pH Analysis and Related Measurements; Approved Guideline – Second Edition. CLSI document C46-A2. Clinical and Laboratory Standards Institute, 950 West Valley Road, Suite 2500, Wayne, Pennsylvania 19087 USA, 2009. 6. Toffaletti J. Effect of small air bubbles on changes in blood pO_2 and blood gas parameters: calculated vs. measured effects. www.acutecaretesting.org Juli 2012. 7. Higgins C. The use of heparin in preparing samples for blood gas analysis. www.acutecaretesting.org April 2007. 8. Siggaard-Andersen O et al. The concentration of free calcium ions in the blood plasma „ionized calcium“. AS-79. Radiometer Medical ApS publication 1980: 163-90. Code no. 918-193. 9. Benoit M et al. Evaluation and advantages of an automatic magnetic mixing of syringes integrated to a whole blood gas analyzer. Scand J Clin Lab Invest 2009; 69(5): 628-32. 10. Grenache D et al. Integrated and automatic mixing of whole blood; an evaluation of a novel blood gas analyzer. Clin Chem Acta 2007; 375: 153-57.

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